

# What our Patients say...

"It's been thrilling from the start to wear *the CRANE* because I never doubted it would work".


**E.D. Wilson, M.D.**

**Discoverer of Wilson's Thyroid Syndrome**

[www.wilsonssyndrome.com](http://www.wilsonssyndrome.com)

"After years of debilitating cranial compression caused by bite equilibration, using *the CRANE* releases the compression in my head and has allowed me to function at a higher level. This device has literally saved my life and brought me from the depths of personal darkness and despair".

Sylvia Kopse  
Ontario



"During the first month of wearing *the CRANE*, a deep posterior pain in my cranium that had plagued me for years vanished. In my disfunction, I became dependent on weekly chiropractic adjustments just to survive. My need for chiropractic adjustments has disappeared since using this device."

Mary K. Fox  
Texas

The Crane N.C.R.P.  
2 Rutledge Circle  
Fletcher, NC 28732

# *the* CRANE

PATENT PENDING



## ELIMINATE Mid-Face / Maxillary Deficiency



The **Only**  
**Non-Cranial**  
Reverse-Pull  
Device

# WHAT IS *the CRANE*?

PATENT PENDING

*the CRANE* is an extra-oral device that stabilizes the head on a cervical collar with a boom for elastics to achieve true forward pull of the maxilla and facial bones.

The weight of the head resting on the cervical collar provides the reactive force to the desired action of the maxillae being pulled forward.

Unlike facemasks, this device revolutionizes forward pull by NOT pushing on the cranial bones, mandible and facial bones to 'achieve results'.

*the CRANE* delivers precise forces to quickly achieve amazing results that were previously unattainable.

*the CRANE* is more comfortable than traditional methods of forward movement because it is not using the mandible (considered an anatomical joint) and facial bones (also subject to being pushed in) to push against to "move" the maxilla forward. This device begs the question do the current models of reverse pull devices really pull the maxilla out OR are the mandible and facial bones being pushed in, thereby creating the illusion of forward movement of the maxilla?

# HOW TO WEAR *the CRANE*

PATENT PENDING

Place *the CRANE* as you would a cervical collar. Adjust the pressure around the neck for the comfort of the patient. The pressure should only be as much as needed to support the head, not to completely immobilize the neck.

Use an orthodontic appliance that will provide the proper support for the elastic force and forward movement.

Call for appliance suggestions.

Sleep with *the CRANE* on !!! The position of the boom and the neck is secured by the collar, so it is safe to sleep with.

## Sizing Guide

.....go to  
[www.thecranencrp.com](http://www.thecranencrp.com)

The grooves on the T-bar of the boom secure elastics for proper direction of forces.



# ADVANTAGES OF *the CRANE*

PATENT PENDING

- Comfortable to wear
- Increased patient compliance
- Forces are NON-CRANIAL

To Order Call

828.243.4131

**\$425.00**

**Includes: Instructions  
Carry Case  
Elastics and Adjusting  
Wrenches**

**Shipping will be an additional cost**

The Crane N.C.R.P., LLC  
2 Rutledge Circle, Fletcher, NC 28732  
828.243.4131  
info@thecranencrp.com  
www.thecranencrp.com